

Why Suffer Pain?

introducing

MULTI PURPOSE MASSAGER

(by MATOBA, JAPAN) since 1962



SELF HELP THERAPY

- Soothing
- Relaxing
- Comforting

Relieve symptoms of stress, ease muscle tension, feel refreshed & revitalised. As used by many professional therapists



Pictures for illustration purposes only

OSKA™ PULSE

PEMF Therapy

PULSE PAIN AWAY

So Easy To Use...



Made in USA

DRUG FREE PAIN & INJURY MANAGEMENT

- A family medical tool suitable for all ages.
- Kick starts the repair processes.
- Accelerates natural healing.

The most portable PEMF Nerve and Muscle stimulator in the world.

More versatile than TENS, Ultrasound and EMS.

just apply to the area of discomfort

(can also be applied through clothing and medical dressings)

Therapy on the go!

From **FREEDOM HEALTH PRODUCTS**

GEOFF HIGGS

www.freedomhealthaustralia.com
3 East Street Maylands WA 6051

GEOFF HIGGS

Phone NOW for your FREE HOME TRIAL!

(08) 9272 5220

Mobile:

0427 11 22 12

Email: higgsy.9876@hotmail.com

RELIEF BEYOND BELIEF

The high-powered multi-directional vibration of the **Matoba** Massager simulates massage by hand. It encourages blood circulation, accelerates metabolism and relaxes muscular tension, a daily tonic massage for about 10 minutes energizes and invigorates the whole body and improves general well-being.

MENTAL TENSION

Massage neck and shoulders with a firm stroking action.

HEADACHE

Massage nape of the neck and lightly on temples.

CONGESTION

Massage deeply in chest area.

FLATULENCE

Massage firmly over abdomen and between shoulder blades.

INSOMNIA

A general tonic massage of the arms, legs, chest, back and abdomen as well as on the feet.

CONSTIPATION

Massage abdominal region for about 3 minutes.

MUSCULAR RHEUMATISM, LUMBAGO, SCIATICA, FIBROSITIS, NEURITIS

Massage knotted muscles and surrounding areas to break down congestion. Use a circular motion for 6-8 minutes

PERIOD PAIN

Massage over lower portion of the back.

MUSCULAR TENSION

Lightly stroke over the painful area for 5 minutes.

ARTHRITIC RHEUMATISM



Massage ligaments and surrounding muscles.

SPRAINS AND CRAMPS

Massage for 3 minutes. Repeat 2-3 times daily

ACHING FEET

Massage the ball and the heel of the foot as well as the ankle for 3 minutes each foot.



Common Vibrators
* IMPACT *



Matoba Action
* HORIZONTAL + VERTICAL + CIRCULAR *

www.freedomhealthaustralia.com
Email: higgsy.9876@hotmail.com

GEOFF HIGGS
FREEDOM HEALTH
PRODUCTS
3 EAST ST., MAYLANDS 6051
PH: 9272 5220
OR 9271 0797
Mob. 0427 11 22 12

Unlike most vibrating units available which employ a thumping/pounding action, **Matoba** Action Massage creates waves of circulating motion that penetrate smoothly and deeply. Used according to directions **Matoba** Action may assist in the maintenance of peripheral circulation